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Breath is Life
we can help you

MPVTRUMA

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SonoDrop 2 Ultrasonic Nebulization
with closed Isapak[®] sterile water system

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Isapak®-Combi-Kit for long-time ultrasonic nebulization of up to 6 days (Part-No. M 50050-07200)

For prevention, alleviation and therapy of respiratory tract disorders

Hygienic safety for the patient is the prime objective

Conventional ultrasonic nebulizers are frequently contaminated and the risk of infection is always imminent.

A **closed sterile water system** such as the *Isapak®* system constitutes a sure and safe solution of the problem. The risk of infection through water-germs is reduced decisively. The aerosol that reaches the patient is produced in a sterile *Isapak®* aerosol chamber. The *Isapak®* aerosol chamber is automatically refilled through an *Isapak®* adapter, also called „chicken-feeder“, which is likewise supplied in a sterile packing. Water is supplied to the system from a 1000 ml sterile water bottle, which has a practical, re-closeable screw cap.

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Isapak®-USN-Kit for short-term ultrasonic nebulization (Part-No. M 55200-00)

Gain a technical lead with SonoDrop 2

Long-time nebulization

Thanks to its high-quality quartz, SonoDrop 2 is particularly suitable for long-time nebulization, e.g. with bedridden patients breathing through the mouth during the night (uninterrupted application up to 8 hours).

Removable quartz

The quartz can be simply removed from the ultrasonic cup with the aid of an assembly device, e.g. in order to also autoclave the ultrasonic cup.

Hygienically safe

Autoclavable ultrasonic cup (without quartz) for still more hygienic safety, e.g. for MRSA-patients.

Simple handling

SonoDrop 2 starts up immediately at the „touch of a button“. The closed *Isapak®* sterile water system is assembled easily in 3 simple steps. This plays an important role with regard to time-management in patient-care.

Integrated heating

A heater is part of the standard features of every SonoDrop 2 basic unit. Simply plug the connection plug of the silicone-heating tube in, and you have all the warm aerosol-mist you need, e.g. for warm air humidification for tracheotomy patients.

Water-refill warning signal

To protect the quartz from overheating and to ensure better therapeutic safety, an acoustic warning signal coupled with a visual signal (red light-emitting diode) will indicate an insufficient water-supply.

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Table-mount, mobile stand or wall-rack mount – the right solution for every field of application!

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Table-mount

The SonoDrop 2 Ultrasonic Nebulizer in table-mounted version is particularly suitable for use at home or in doctor's practices. For fast moisturization of the respiratory tract to alleviate acute discomfort (e.g. in case of pseudocroup - laryngitis) at home and primarily in the doctor's practice, therapy with SonoDrop 2 is a very effective natural method of treatment.

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Wall-rack mount

SonoDrop 2 can be universally utilized, e.g. in clinical operation mounted on wall-racks. With the special wall-rack mount, SonoDrop 2 can be ideally adapted.

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The application must always be prescribed by a physician.

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Mobile stand

The standard mobile stand is suitable for demanding service in the clinic, in the nursing home or in the doctor's practice. The standard mobile stand is made of sturdy stainless steel, and it is equipped with durable five-caster pedestal with locking brakes.

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Technical Data:

<i>Product Class:</i>	<i>II a according to MPG (Act on Medical Devices)</i>
<i>Type:</i>	<i>SonoDrop Ultrasonic Nebulizer</i>
<i>Measurements basic unit:</i>	<i>approx. height 140 mm, width 220 mm, depth 300 mm</i>
<i>Weight basic unit / total:</i>	<i>3.5 kg/ 6 kg</i>
<i>Nebulizing performance:</i>	<i>up to 2.5 ml/min</i>
<i>Flow rate:</i>	<i>min. 10.2 l/min ± 0.4 l/min max. 18.7 l/min ± 0.4 l/min</i>
<i>Size of aerosol particles:</i>	<i>1 - 12 µm, 4.5 µm MMD (Median Mass Diameter, measured with Malvern Laser Particle Sizer at 23°C) 50 % smaller than 5 µm 23.5 % smaller than 3 µm 17.0 % smaller than 2 µm 12.0 % smaller than 1 µm</i>
<i>Measurements standard mobile stand:</i>	<i>weight: 13.0 kg height: 1.35 m base surface: ø 0.60 cm</i>
<i>Nominal Voltage</i>	<i>230 V AC 50 Hz</i>
<i>Mains fuse</i>	<i>2 x 1 A, delayed action (slo-blo)</i>
<i>Nominal current without / with heating</i>	<i>0.27 A / 0.4 A</i>
<i>Power consumption without / with heating</i>	<i>62 VA / 92 VA</i>
<i>HF-Frequency</i>	<i>1.65 MHz</i>
<i>Sound pressure level:</i>	<i>approx. 26 dbA (silent)</i>
<i>Operating Mode</i>	<i>continuous operation</i>
<i>Protection Class</i>	<i>I, type B</i>
<i>Aerosol temperature</i>	<i>up to 37 °C, if heating tube is used</i>
<i>Interference suppression</i>	<i>EN 60601-1-2</i>
<i>Operating conditions:</i>	<i>+15 °C to +40 °C 30 - 70 % relative humidity The tube heating may not be used at environmental temperatures of more than 30 °C</i>
<i>Caution:</i>	<i>In accordance with the EU Directive on Medical Products (93/42 EWG) and the Act on Medical Devices(MPG)</i>
CE 0123	
<i>Guarantee:</i>	<i>12 months as of date of purchase</i>
<i>Manufacturer:</i>	<i>MPV TRUMA GmbH Wernher-von-Braun-Str. 12 D 85640 Putzbrunn</i>

Return device to manufacturer for disposal. Subject to technical amendments!

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SonoDrop 2 plus ideal for operation in stationary facilities

*SonoDrop 2 is also available in a „plus“ version, especially suited for use in the hospital, the doctor's practice or in geriatric care. Besides the proven, reliable functions, **SonoDrop 2 plus** features a „plus“ of parameters, which are of importance for daily operation in hospitals, in geriatric care centers or in doctor's practices (see separate brochure).*

- ▶ **Timer**
Various intervals can be pre-set with an automatic timer.
- ▶ **Flow adjustment**
- ▶ **Aerosol temperature**
- ▶ **Nebulizing performance / aerosol density**

With the digital display, all parameters can be set and checked safely and with ease. An integrated running time meter provides for complete control regarding Compliance.

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A well-arranged display ensures simple setting and adjustment of all parameters.

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Handling of the Isapak[®] system – captivatingly simple

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Isapak[®] USN-Kombi-Kit

The economic Isapak[®] set for long-time application, consisting of:
6 x Isapak[®] sterile water, 1000 ml
1 x Isapak[®] aerosol chamber
1 x Bacterial filter
1 x Isapak[®] adapter

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Step 1

Remove the Isapak[®] aerosol chamber from the sterile packing and position it on the ultrasonic cup. Make a note of the date the pack was opened on the Isapak[®] aerosol chamber.

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Step 2

Place the Isapak[®] sterile water bottle on the table and open the screw-cap. Remove the Isapak[®] adapter from the sterile packing, and screw it onto the sterile water bottle. Note: The Isapak[®] adapter is autoclavable (up to 121°C) and it can be used repeatedly. Record the date the pack was opened on the Isapak[®] sterile water bottle.

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Step 3

Introduce the Isapak[®] sterile water bottle with the adapter from above through the bottle holder into the opening of the Isapak[®] aerosol chamber.

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
Innovative technology for even more reliability

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Lufteinlassfilter / *Air inlet filter*

Bakterienfilter / *Bacterial filter*

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Integrated heating

Every basic unit is generally equipped with a heating unit. An ultra-light silicone heating tube provides for an aerosol temperature of approx. 37 °C.

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New ultrasonic cup concept

Autoclavable ultrasonic cup

The SonoDrop 2 – ultrasonic cup can be taken apart completely. I.e., the quartz can be removed from the ultrasonic cup using an assembly device in order to be able to autoclave the ultrasonic cup and the quartz receptacle at 134 °C. Please do not autoclave the quartz!

Automatic water detection

If there is no water in the ultrasonic cup, the device will not start, therefore protecting the quartz from overheating. This is indicated visually by a yellow light-emitting diode at the front of the device. An acoustic signal additionally signifies this condition.

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Clean inhalation air

*An air inlet filter located on the bottom side of the SonoDrop 2 basic unit provides for a preliminary filtration of the dust particles contained in the intake-air. A bacteria- and virus-filter BARR-VENT-S – with a bacteria retention capacity of 99.988% and with a reduction factor of 9.3×10^3 for *Pseudomonas aeruginosa* – provides for a hygienically safe inhalation therapy.*